# Body Screen TX, weiß

Art.-Nr.: 75001


**HAMMER**

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PACKAGE LIST

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Body Screen TX, weiß

**4 AA Batterien**  
4 AA Batteries  
4 piles AA

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# Table of content

01	Important Safety Guidelines	P. 10
02	Specifications	P. 11
03	Getting to Know Your Body Screen TX	P. 12
04	Display	P. 13
05	App Installation	P. 13
06	Using Your Body Screen TX	P. 14
07	Care & Maintenance	P. 16
08	Warranty	P. 16

# Important Safety Guidelines

When using electrical appliance, basic safety precautions should always be followed, including the following:

Read all instructions carefully before using the Body Screen TX. Retain this user manual for future reference.

## WARNINGS

- Individuals with medical implant devices such as pacemakers, or essential support devices such as patient monitoring systems, must not use this equipment. Safe, low-level currents will flow through the body during the analysis, which may cause malfunctioning of the device or endanger lives.
- Bioelectrical impedance analysis (BIA) uses safe low-current values that are not harmful to the body. However, we advise pregnant women not to undergo an analysis.
- The surface of the Body Screen TX can become slippery in wet places. Please watch your step and do not slip.

## CAUTION

- Stand upright for about 5 minutes before the analysis. If you carry out the analysis immediately after lying down or sitting for a long time, the analysis results may change slightly. This is because body water tends to move to the lower body as soon as the person stands or gets up.
- Do not eat before the analysis. In cases where the examinee has already eaten, the analysis should be put off for at least two hours after the meal. This is because food mass is included in the examinee's weight and thus, may result in measurement errors.
- Use the bathroom before the analysis. Waste is not included in the body's compositional elements, but the volume of urine and excrement is included in the weight measurement affecting accuracy of the analysis results.
- Do not exercise before the analysis. Strenuous exercise or sharp movements can cause temporary changes in body composition. Even light exercise can change your body composition temporarily.
- Take the analysis in the morning, if possible. Body water tends to gravitate towards the lower body throughout the day, affecting accuracy of the analysis results.
- Take the analysis at the same time and same condition in each day to reduce test variables, if possible, to achieve a coherent analysis record.
- Thoroughly wipe the palms and soles with wet tissue/towel before the analysis. The analysis may be difficult if the examinee's palms and soles are too dry or if the examinee has too many calluses.
- Put the HAMMER Body Screen TX on a leveled, hard and non-vibrating surface. Put the scale on an uneven surface may cause the examinee to fall down. Analysis results may also be inaccurate.
- Avoid contact with the examinee during the analysis. Contact may lead to interference affecting analysis results.


# Specifications

Bioelectrical Impedance Analysis (BIA)	15 Impedance Measurements by Using 3 Different Frequencies ( 5kHz, 50kHz,250kHz )at Each of 5 Segment ( Right Arm, Left Arm, Trunk, Right Leg, and Left Leg )
Electrode	Tetrapolar 8-Point Tactile Electrodes
Measurement Method	Direct Segmental Multi-frequency Bioelectrical Impedance Analysis (DSM-BIA), Simultaneous Multi-frequency Impedance Measurement (SMFIM)
43 Body Analyses	Total weight, ideal weight, difference from ideal weight, total lean mass, basal metabolic rate, BMI, body fat percentage, total body water, body water percentage, total protein, visceral fat, total minerals, bone minerals, total muscle mass, difference from ideal muscle mass, skeletal muscle mass, segmental muscle mass (left arm, right arm, trunk, left leg, right leg), body fat mass, Difference from ideal fat mass, segmental fat (left arm, right arm, trunk, left leg, right leg), segmental lean mass, waist-to-hip ratio, body build, edema evaluation, circumference of body parts (neck, chest, waist, hip, right arm, left arm, right thigh, left thigh), body balance (upper balance, lower balance, upper-lower balance), total score.
Database	The firmware can save up to 10 results. The App can save up unlimited results.
Connection	Bluetooth 5.0
Battery Type	4 AA Batteries
Battery Life	Around 3 months
Display Type	Hidden LED Screen
Dimension	31.5 x 31.5 x 3.3 cm
Duration of analysis	About 20 seconds
Weight	2.35 kg
Operation Environment	10 ~ 40°C
Storage Environment	-10 ~ 60°C
User weight	5 ~ 200kg
User age	13 ~ 99 years

# Get To Know your HAMMER Body Screen TX


HAMMER Body Screen TX detects and analyzes body composition of different body parts by BIA (Bioelectrical Impedance Analysis). It uses 8-point tactile electrodes, multiple currents at low, medium and high frequencies, direct measurement of each of the five body parts (directly measure each of the body's five cylinders left arm, right arm, torso, left leg, and right leg) and the basic properties of the human body to deliver quick, accurate and detailed results.

# Product Diagram

# Display

LED EFFECT	OFF	ON	Flash
	Not connected	Connected	Not connected
	Not measuring	Measuring	/
	Not connected	Connected	/
	Not connected	Connected	/
	Not measuring	Measuring	/

# App Installation

Search "BodyPedia" in the Android Application Market or App Store to download and install the App.

App downloaded via other channels may be faulty. Before pairing with your equipment, make sure your equipment meets the following requirements:

- iPhone 5S and above with iOS 10.0 and above;
- Android 5.0 and above with Bluetooth 4.0 and above

Connect to the Body Screen TX via Bluetooth

Bind the Body Screen TX to BodyPedia App


Open BodyPedia App and Bluetooth of the phone, follow the instruction of the App to register, bind the Body Screen TX to the App and finish the first measurement.

Bluetooth connection


Step on the scale with one foot to activate the Body Screen TX, and take your foot off the Body Screen TX when the scale screen lights up. Open BodyPedia App, and the Body Screen TX will connect to the App automatically if the scale has been bound to the App.

# Using Your Body Screen TX

1. Pull out the tab on the underside of the Body Screen TX. Place the Body Screen TX on the floor and switch it on by briefly putting weight on it with one foot. The screen will now display 0.0 kg/lb.


2. Now stand barefoot on the Body Screen TX. The weight measurement begins..

Stand calmly and centered on the Body Screen TX.


The Body Screen TX must be placed on a level and firm surface, e.g. on parquet flooring or a tile.

3. As soon as the weight stops flashing, the weight measurement is complete.


The number stops flashing.

4. The Body Screen TX is now ready for body analysis. Lean towards the handles and pull them out with one click.

5. The body test begins.

The examinee must maintain proper posture to have accurate analysis results.

6. When the test is completed, connect the Body Screen TX to your App to have the results.

Feet			
Legs			
Arms			
Hands			

## Care & Maintenance

Do not extend the handles beyond its limitation.

Do not place any objects on the footplate.

Do not apply excessive force on the equipment.

Do not allow any liquid substances to contact the equipment directly. Keep food and drinks away from the equipment. Substances getting inside the equipment can cause critical damage to the electronic components.

Remove the batteries when they become weak or when the body analyzer is not in use for extended periods. This will prevent damage to the device from battery leakage at such times

Use a lint-free cloth to gently wipe the external surface of the equipment about once every week

Be careful not to scratch the LED screen.

HAMMER Body Screen TX does not need regular maintenance. If some problems occur while operating the device, get in touch with our customer service team.

We do not take the responsibility about problems caused by any arbitrary repairs.

## Warranty

We provide warranty for our devices according to the following conditions:

- 1) We will repair free of charge in accordance with the following conditions: Damage or defects to the device which are demonstrably due to a manufacturing defect, if they are reported to us immediately upon discovery and within 24 months after delivery to the end user. The warranty does not cover easily breakable parts, such as glass or plastic. A warranty obligation is not triggered: by minor deviations from the target condition, which are insignificant for the value and usability of the device, by damage from chemical and electrochemical effects, from water and generally by damage due to force majeure.
- 2) The warranty service is provided in such a way that defective parts are repaired free of charge or replaced by faultless parts at our discretion. The costs for material and working time will be borne by us. Repairs at the place of installation cannot be demanded. The proof of purchase with the date of purchase and/or delivery must be presented. Replaced parts become our property.
- 3) The warranty claim expires if repairs or interventions are carried out by persons who are not authorized by us to do so or if our devices are equipped with supplementary or accessory parts which are not matched to our devices, furthermore if the device is damaged or destroyed by the influence of force majeure or by environmental influences, in the case of damage which has occurred due to improper handling - in particular non-observance of the operating instructions - or maintenance or if the device shows mechanical damage of any kind. The customer service can authorize you after telephone consultation to repair or exchange parts, which will be sent to you. In this case, the warranty claim will of course not expire.
- 4) Warranty services shall neither cause an extension of the warranty period nor shall they set a new warranty period in motion.
- 5) Further or other claims, in particular those for compensation for damage occurring outside the device, are excluded - insofar as liability is not mandatorily ordered by law.
- 6) Our warranty conditions, which include the prerequisites and scope of our warranty services, do not affect the contractual warranty obligations of the seller.
- 7) Parts subject to wear and tear, which are identified as such in the parts list, are not subject to the warranty conditions.
- 8) The warranty claim expires in case of improper use, especially in fitness studios, rehabilitation facilities and hotels.

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01	Wichtige Sicherheitsrichtlinien	P. 01
02	Spezifikationen	P. 02
03	Lerne deine HAMMER Body Screen TX kennen	P. 03
04	Display	P. 04
05	App-Installation	P. 04
06	Verwendung der Body Screen TX	P. 05
07	Pflege & Wartung	P. 07
08	Garantiebedingungen	P. 07
09	Service-Kontakt	P. 08

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Hammer Sport AG